

BHMT Rabbit mAb

Catalog No: #56665

Package Size: #56665-1 50ul #56665-2 100ul


Description

Product Name	BHMT Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human
Specificity	BHMT Antibody detects endogenous levels of total BHMT
Immunogen Description	A synthesized peptide derived from human BHMT
Conjugates	Unconjugated
Other Names	BHMT; BHMT1; HEL S 61p;
Accession No.	Uniprot:Q93088
Calculated MW	45kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200IP:1:50

Images



Western blot analysis of BHMT expression in Human fetal kidney lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.