

Recombinant IL6 Rabbit mAb(G12)

Catalog No: #58001

Package Size: #58001-1 50ul #58001-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Recombinant IL6 Rabbit mAb(G12)
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	G12
Purification	Affinity purification
Applications	WB,IHC,ELISA
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total IL6 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human IL6
Conjugates	Unconjugated
Target Name	IL6
Other Names	IL6;BSF-2;BSF2;CDF;HGF;HSF;IFN-beta-2;IFNB2;IL-6
Accession No.	Uniprot:P05231GeneID:3569
SDS-PAGE MW	26KDa
Concentration	0.7mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details


WB□1:500 - 1:2000

IHC□1:50 - 1:200


Images



Immunohistochemical analysis of paraffin-embedded human lung cancer tissue



Immunohistochemical analysis of paraffin-embedded human lung cancer tissue



Human IL-6 ELISA Standard Curve/ Detection Antibody
 oO G12oO oO 58001

Note: This product is for in vitro research use only and is not intended for use in humans or animals.