

CCR2/CKR2 Rabbit mAb

Catalog No: #58558

Package Size: #58558-1 50ul #58558-2 100ul


Description

Product Name	CCR2/CKR2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	IHC FC IP
Species Reactivity	Human Mouse
Specificity	CCR2/CKR2 Antibody detects endogenous levels of total CCR2/CKR2
Immunogen Description	A synthesized peptide derived from human CCR2/CKR2
Other Names	C-C chemokine receptor type 2; C-C CKR-2; CC chemokine receptor type 2; CC CKR 2; CC-CKR-2; CCCKR2; CCR 2; CCR-2;; CCR1L;; CCR2;
Accession No.	Uniprot:P41597
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

IHC 1:50~1:200 IP 1:50 FC1:100

Images



Immunohistochemical analysis of paraffin-embedded human tonsil, using CCR2/CKR2 Antibody.

Background

Receptor for the MCP-1, MCP-3 and MCP-4 chemokines. Transduces a signal by increasing the intracellular calcium ions level. Alternative coreceptor with CD4 for HIV-1 infection.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.