

Tissue Factor Rabbit mAb

Catalog No: #59064

Package Size: #59064-1 50ul #59064-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Tissue Factor Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	Tissue Factor Antibody detects endogenous levels of Tissue Factor
Immunogen Description	A synthesized peptide derived from human Tissue Factor
Other Names	CD142 antigen; Coagulation factor III (thromboplastin tissue factor); F3; TF; TFA; Thromboplastin; Tissue factor;
Accession No.	Uniprot:P13726
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200

Images



Western blot analysis of Tissue Factor expression in SHSY5Y cell lysate.

Background

Initiates blood coagulation by forming a complex with circulating factor VII or VIIa. The [TF:VIIa] complex activates factors IX or X by specific limited proteolysis. TF plays a role in normal hemostasis by initiating the cell-surface assembly and propagation of the coagulation protease cascade.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.