

Hsp75 Rabbit mAb

Catalog No: #59176

Package Size: #59176-1 50ul #59176-2 100ul


Description

Product Name	Hsp75 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human
Specificity	Hsp75 Antibody detects endogenous levels of total Hsp75
Immunogen Description	A synthesized peptide derived from human Hsp75
Other Names	HSP75; HSP 75; HSP90L; TRAP-1;
Accession No.	Uniprot:Q12931
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IP 1:50

Images



Western blot analysis of Hsp75 expression in K562 cell lysate.

Background

TNF receptor-associated protein 1 (TRAP1), also known as HSP75, is a mitochondrial chaperone and ATPase that was originally identified as a protein that interacts with the TNF receptor.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.