

Transferrin Receptor (CD71) Rabbit mAb

Catalog No: #59241

Package Size: #59241-1 50ul #59241-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Transferrin Receptor (CD71) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP FC
Species Reactivity	Human Mouse Rat
Specificity	TFR Antibody detects endogenous levels of total TFR
Immunogen Description	A synthesized peptide derived from human TFR
Other Names	CD71 antigen; T9; TFRC; TR; TfR1; Transferrin receptor protein 1; Trfr; p90;
Accession No.	Uniprot:P02786
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB 1:500~1:2000 IHC 1:50~1:200 IP 1:50 FC 1:50

Images



Western blot analysis of Transferrin Receptor (CD71) expression in HeLa cell lysate.



Immunohistochemical analysis of paraffin-embedded human placenta, using Transferrin Receptor (CD71) Antibody.

Background

Binds the telomeric double-stranded TTAGGG repeat and negatively regulates telomere length. Involved in the regulation of the mitotic spindle. Component of the shelterin complex (telosome) that is involved in the regulation of telomere length and protection.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.