

RAGE Rabbit mAb

Catalog No: #59253

Package Size: #59253-1 50ul #59253-2 100ul


Description

Product Name	RAGE Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	RAGE Antibody detects endogenous levels of total RAGE
Immunogen Description	A synthesized peptide derived from human RAGE
Other Names	Advanced glycosylation end product-specific receptor; Ager;
Accession No.	Uniprot:Q15109
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of RAGE expression in mouse lung lysate.

Background

The receptor for advanced glycation end products (RAGE) is member of the immunoglobulin (Ig) superfamily. Mediates interactions of advanced glycosylation end products (AGE). Binding of AGEs to RAGE results in the induction of cellular oxidant stress and activation of the transcription factor NFkB. Evidence suggests that the induction of oxidant stress results in the activation of an intracellular cascade involving p21 ras and MAP kinase, which leads to activation of transcription.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.