

Calbindin Rabbit mAb

Catalog No: #59325

Package Size: #59325-1 50ul #59325-2 100ul


Description

Product Name	Calbindin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	Calbindin Antibody detects endogenous levels of total Calbindin
Immunogen Description	A synthesized peptide derived from human Calbindin
Other Names	avian-type; CAB27; CALB1; Calbindin 1 28kDa; Calbindin; D28K; RTVL H protein;
Accession No.	Uniprot:P05937
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200

Images



Western blot analysis of Calbindin expression in human fetal brain lysate.

Background

Buffers cytosolic calcium. May stimulate a membrane Ca(2+)-ATPase and a 3',5'-cyclic nucleotide phosphodiesterase.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.