

CD32 Rabbit mAb

Catalog No: #59339

Package Size: #59339-1 50ul #59339-2 100ul


Description

Product Name	CD32 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	25C8
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	CD32 Antibody detects endogenous levels of total CD32
Immunogen Description	A synthesized peptide derived from human CD32
Target Name	CD32
Other Names	CD32; CD32A; CD32B; CD32C; FCG2; FcGR; FCGR2; FCGR2A; FCGR2B; FCGR2C; Fcgr3a; FCRII; IGFR2;
Accession No.	Uniprot:P12318
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Images



All Lanes:CD32 Antibody at 1/ 1000 dilution.Lane1:K562 cell lysate.Proteins at 20ug per laneSecondary: Goat Anti-Rabbit IgG(HRP) at 1/20000 dilutionPredicted band size : 43kDaObserved band size: 43kDa

Background

Binds to the Fc region of immunoglobulins gamma. Low affinity receptor. By binding to IgG it initiates cellular responses against pathogens and soluble antigens.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.