

NUR77 Rabbit mAb

Catalog No: #59458

Package Size: #59458-1 50ul #59458-2 100ul


Description

Product Name	NUR77 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	NUR77 Antibody detects endogenous levels of total NUR77
Immunogen Description	A synthesized peptide derived from human NUR77
Other Names	Nur77; Orphan nuclear receptor HMR; Orphan nuclear receptor TR3; ST-59; Testicular receptor 3; NR4A1; GFRP1; HMR; NAK1;
Accession No.	Uniprot:P22736
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of NUR77 expression in HepG2 cell lysate.

Background

Orphan nuclear receptor. May act concomitantly with NURR1 in regulating the expression of delayed-early genes during liver regeneration. Binds the NGFI-B response element (NBRE) 5'-AAAAGGTCA-3' (By similarity). May inhibit NF-kappa-B transactivation of IL2.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.