

## BST2 Rabbit mAb

Catalog No: #59498

Package Size: #59498-1 50ul #59498-2 100ul


## Description

|                       |                                                                                                    |
|-----------------------|----------------------------------------------------------------------------------------------------|
| Product Name          | BST2 Rabbit mAb                                                                                    |
| Host Species          | Rabbit                                                                                             |
| Clonality             | Monoclonal                                                                                         |
| Isotype               | Rabbit IgG                                                                                         |
| Purification          | Affinity-chromatography                                                                            |
| Applications          | WB IHC FC                                                                                          |
| Species Reactivity    | Human                                                                                              |
| Specificity           | BST2 Antibody detects endogenous levels of total BST2                                              |
| Immunogen Description | A synthesized peptide derived from human BST2                                                      |
| Other Names           | Bone marrow stromal antigen 2; BST2; CD317; HM1.24 antigen; NPC A 7; Tetherin;                     |
| Accession No.         | Uniprot:Q10589                                                                                     |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage               | Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.                     |

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200 FC 1:50

## Images



Western blot analysis of BST2 expression in HeLa cell lysate.

## Background

May be involved in the sorting of secreted proteins (By similarity). May be involved in pre-B-cell growth. Antiretroviral defense protein, that blocks release of retrovirus from the cell surface.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.