

TCTP Rabbit mAb

Catalog No: #59510

Package Size: #59510-1 50ul #59510-2 100ul


Description

Product Name	TCTP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	TCTP Antibody detects endogenous levels of total TCTP
Immunogen Description	A synthesized peptide derived from human TCTP
Other Names	Fortilin; Histamine releasing factor; HRF; Immunoglobulin E-dependent; p02; p23; TCTP; tpt1;
Accession No.	Uniprot:P13693
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:5000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50 FC 1:100

Images



Western blot analysis of TCTP Receptor expression in HeLa cell lysate.

Background

Involved in calcium binding and microtubule stabilization.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.