

SNAP23 Rabbit mAb

Catalog No: #59523

Package Size: #59523-1 50ul #59523-2 100ul


Description

Product Name	SNAP23 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	SNAP23 Antibody detects endogenous levels of total SNAP23
Immunogen Description	A synthesized peptide derived from human SNAP23
Other Names	HsT17016; LS-B8340; SNAP23A; SNAP23B; Synaptosomal associated protein 23; Vesicle membrane fusion protein SNAP23;
Accession No.	Uniprot:O00161
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of SNAP23 expression in HeLa cell lysate.

Background

Essential component of the high affinity receptor for the general membrane fusion machinery and an important regulator of transport vesicle docking and fusion.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.