

Reelin Rabbit mAb

Catalog No: #59529

Package Size: #59529-1 50ul #59529-2 100ul


Description

Product Name	Reelin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Rat
Specificity	Reelin Antibody detects endogenous levels of total Reelin
Immunogen Description	A synthesized peptide derived from human Reelin
Other Names	LIS2; PRO1598; Reeler; Reelin; RELN; RL;
Accession No.	Uniprot:P78509
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of Reelin expression in HepG2 cell lysate.

Background

Extracellular matrix serine protease that plays a role in layering of neurons in the cerebral cortex and cerebellum. Regulates microtubule function in neurons and neuronal migration. Affects migration of sympathetic preganglionic neurons in the spinal cord, where it seems to act as a barrier to neuronal migration.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.