

Fragilis Rabbit mAb

Catalog No: #59574

Package Size: #59574-1 50ul #59574-2 100ul


Description

Product Name	Fragilis Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human
Specificity	Fragilis Antibody detects endogenous levels of total Fragilis
Immunogen Description	A synthesized peptide derived from human Fragilis
Other Names	Fragilis; IFITM3; Interferon inducible; IP15;
Accession No.	Uniprot:Q01628
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of Fragilis expression in HeLa cell lysate.

Background

IFN-induced antiviral protein that mediates cellular innate immunity to at least three major human pathogens, namely influenza A H1N1 virus, West Nile virus (WNV), and dengue virus (WNV), by inhibiting the early step(s) of replication.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.