

USP28 Rabbit mAb

Catalog No: #59589

Package Size: #59589-1 50ul #59589-2 100ul


Description

Product Name	USP28 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human
Specificity	USP28 Antibody detects endogenous levels of total USP28
Immunogen Description	A synthesized peptide derived from human USP28
Other Names	KIAA1515; Ubiquitin carboxyl terminal hydrolase 28; Ubiquitin thioesterase 28; UBP28;
Accession No.	Uniprot:Q96RU2
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of USP28 expression in A431 cell lysate.

Background

Deubiquitinase involved in DNA damage response checkpoint and MYC proto-oncogene stability. Involved in DNA damage induced apoptosis by specifically deubiquitinating proteins of the DNA damage pathway such as CLSPN.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.