

STIP1 Rabbit mAb

Catalog No: #59630

Package Size: #59630-1 50ul #59630-2 100ul


Description

Product Name	STIP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	STIP1 Antibody detects endogenous levels of total STIP1
Immunogen Description	A synthesized peptide derived from human STIP1
Other Names	Hop; P60; STI1; STI1L; STIP1; Stress induced phosphoprotein 1;
Accession No.	Uniprot:P31948
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:5000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50 FC 1:80

Images



Western blot analysis of STIP1 expression in Jurkat cell lysate.

Background

Mediates the association of the molecular chaperones HSC70 and HSP90 (HSPCA and HSPCB).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.