

ALIX Rabbit mAb

Catalog No: #59649

Package Size: #59649-1 50ul #59649-2 100ul


Description

Product Name	ALIX Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse Rat
Specificity	ALIX Antibody detects endogenous levels of total ALIX
Immunogen Description	A synthesized peptide derived from human ALIX
Other Names	AIP1; Alix; DRIP4; Hp95; PDCD6IP;
Accession No.	Uniprot:Q8WUM4
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 FC 1:100

Images



Western blot analysis of ALIX expression in K562 cell lysate.

Background

Class E VPS protein involved in concentration and sorting of cargo proteins of the multivesicular body (MVB) for incorporation into intraluminal vesicles (ILVs) that are generated by invagination and scission from the limiting membrane of the endosome.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.