

MMP17 Rabbit mAb

Catalog No: #59660

Package Size: #59660-1 50ul #59660-2 100ul


Description

Product Name	MMP17 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	MMP17 Antibody detects endogenous levels of total MMP17
Immunogen Description	A synthesized peptide derived from human MMP17
Other Names	MMP 17; Mmp17; MT MMP4; MT-MMP 4; MT4 MMP; MT4-MMP; MT4MMP; MTMMP4;
Accession No.	Uniprot:Q9ULZ9
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:1000 IHC 1:50~1:200 ICC/IF 1:50~1:200

Images



Western blot analysis of MMP17 expression in HeLa cell lysate.

Background

Endopeptidase that degrades various components of the extracellular matrix, such as fibrin. May be involved in the activation of membrane-bound precursors of growth factors or inflammatory mediators, such as tumor necrosis factor-alpha. May also be involved in tumoral process.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.