

NR0B1 / DAX1 Rabbit mAb

Catalog No: #59661

Package Size: #59661-1 50ul #59661-2 100ul


Description

Product Name	NR0B1 / DAX1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human
Specificity	NR0B1 / Dax1 Antibody detects endogenous levels of total NR0B1 / Dax1
Immunogen Description	A synthesized peptide derived from human NR0B1 / Dax1
Other Names	AHC; AHX; DSS; GTD; HHG; AHCH; DAX1; DAX-1; NROB1; SRXY2;
Accession No.	Uniprot:P51843
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:1000 IHC 1:100~1:500 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of NR0B1 / Dax1 expression in A431 cell lysate.

Background

Orphan nuclear receptor. Component of a cascade required for the development of the hypothalamic-pituitary-adrenal-gonadal axis. Acts as a coregulatory protein that inhibits the transcriptional activity of other nuclear receptors through heterodimeric interactions.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.