

GIT1 Rabbit mAb

Catalog No: #59671

Package Size: #59671-1 50ul #59671-2 100ul


Description

Product Name	GIT1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse Rat
Specificity	GIT1 Antibody detects endogenous levels of total GIT1
Immunogen Description	A synthesized peptide derived from human GIT1
Other Names	ARF GAP GIT1; CAT 1; CaT1; GIT1;
Accession No.	Uniprot:Q9Y2X7
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 FC 1:100

Images



Western blot analysis of GIT1 expression in 293T cell lysate.

Background

GTPase-activating protein for the ADP ribosylation factor family. May serve as a scaffold to bring together molecules to form signaling modules controlling vesicle trafficking, adhesion and cytoskeletal organization. Increases the speed of cell migration, as well as the size and rate of formation of protrusions, possibly by targeting PAK1 to adhesions and the leading edge of lamellipodia.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.