

Syntrophin alpha 1 Rabbit mAb

Catalog No: #59739

Package Size: #59739-1 50ul #59739-2 100ul


Description

Product Name	Syntrophin alpha 1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse
Specificity	Syntrophin alpha 1 Antibody detects endogenous levels of total Syntrophin alpha 1
Immunogen Description	A synthesized peptide derived from human Syntrophin alpha 1
Other Names	Alpha 1 syntrophin; LQT12; SNT1; Snta1; Syntrophin 1; TACIP1;
Accession No.	Uniprot:Q13424
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:100

Images



Western blot analysis of Syntrophin alpha 1 expression in A549 cell lysate.

Background

Adapter protein that binds to and probably organizes the subcellular localization of a variety of membrane proteins. May link various receptors to the actin cytoskeleton and the extracellular matrix via the dystrophin glycoprotein complex.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.