

Cdx1 Rabbit mAb

Catalog No: #59825

Package Size: #59825-1 50ul #59825-2 100ul


Description

Product Name	Cdx1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human
Specificity	Cdx1 Antibody detects endogenous levels of total Cdx1
Immunogen Description	A synthesized peptide derived from human Cdx1
Other Names	CDX1; Homeobox protein CDX 1;
Accession No.	Uniprot:P47902
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:1000 ICC/IF 1:50~1:200

Images



Western blot analysis of Cdx1 expression in Caco-2 cell lysate.

Background

Could play a role in the terminal differentiation of the intestine.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.