

## ADAMTS13 Rabbit mAb

Catalog No: #59847

Package Size: #59847-1 50ul #59847-2 100ul


## Description

|                       |                                                                                                    |
|-----------------------|----------------------------------------------------------------------------------------------------|
| Product Name          | ADAMTS13 Rabbit mAb                                                                                |
| Host Species          | Rabbit                                                                                             |
| Clonality             | Monoclonal                                                                                         |
| Isotype               | Rabbit IgG                                                                                         |
| Purification          | Affinity-chromatography                                                                            |
| Applications          | WB FC                                                                                              |
| Species Reactivity    | Human                                                                                              |
| Specificity           | ADAMTS13 Antibody detects endogenous levels of total ADAMTS13                                      |
| Immunogen Description | A synthesized peptide derived from human ADAMTS13                                                  |
| Other Names           | ADAMTS13; ADAM-TS; ADAM-TS13; ADAMTS-13; ADAM-TS 13; C9orf8; TTP; VWF-CP; VWFCP;                   |
| Accession No.         | Uniprot:Q76LX8                                                                                     |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage               | Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.                     |

## Application Details

WB 1:500~1:2000 FC 1:100

## Images



Western blot analysis of ADAMTS13 expression in A549 cell lysate.

## Background

Cleaves the vWF multimers in plasma into smaller forms.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.