

GDF15 Rabbit mAb

Catalog No: #59864

Package Size: #59864-1 50ul #59864-2 100ul


Description

Product Name	GDF15 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	GDF15 Antibody detects endogenous levels of total GDF15
Immunogen Description	A synthesized peptide derived from human GDF15
Other Names	GDF-15; growth differentiation factor 15; MIC-1; MIC1; NAG-1; PDF; PLAB; PTGFB;
Accession No.	Uniprot:Q99988
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of GDF15 expression in HepG2 cell lysate.

Background

Highly expressed in placenta, with lower levels in prostate and colon and some expression in kidney.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.