

FAAH1 Rabbit mAb

Catalog No: #59878

Package Size: #59878-1 50ul #59878-2 100ul


Description

Product Name	FAAH1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse
Specificity	FAAH1 Antibody detects endogenous levels of total FAAH1
Immunogen Description	A synthesized peptide derived from human FAAH1
Other Names	FA2H; FAAH; FAAH-1; Oleamide hydrolase 1; Oleamide hydrolase;
Accession No.	Uniprot:O00519
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Images



Western blot analysis of FAAH1 expression in A431 cell lysate.

Background

Degrades bioactive fatty acid amides like oleamide, the endogenous cannabinoid, anandamide and myristic amide to their corresponding acids, thereby serving to terminate the signaling functions of these molecules.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.