

BANF1 Rabbit mAb

Catalog No: #59881

Package Size: #59881-1 50ul #59881-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	BANF1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	BANF1 Antibody detects endogenous levels of total BANF1
Immunogen Description	A synthesized peptide derived from human BANF1
Other Names	BAF; BANF1; BCRG1; BCRP1; NGPS;
Accession No.	Uniprot:O75531
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:100

Images



Western blot analysis of BANF1 expression in HeLa cell lysate.

Background

Plays fundamental roles in nuclear assembly, chromatin organization, gene expression and gonad development. May potentially compress chromatin structure and be involved in membrane recruitment and chromatin decondensation during nuclear assembly.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.