

SNF2H Rabbit mAb

Catalog No: #59882

Package Size: #59882-1 50ul #59882-2 100ul


Description

Product Name	SNF2H Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	SNF2H Antibody detects endogenous levels of total SNF2H
Immunogen Description	A synthesized peptide derived from human SNF2H
Other Names	hISWI; hSNF2H; ISWI; MommeD4; SMARCA5; Snf2h; WCRF135;
Accession No.	Uniprot:O60264
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200

Images



Western blot analysis of SNF2H expression in MCF7 cell lysate.

Background

Helicase that possesses intrinsic ATP-dependent nucleosome-remodeling activity. Complexes containing SMARCA5 are capable of forming ordered nucleosome arrays on chromatin; this may require intact histone H4 tails.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.