

Aconitase 2 Rabbit mAb

Catalog No: #59893

Package Size: #59893-1 50ul #59893-2 100ul


Description

Product Name	Aconitase 2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	Aconitase 2 Antibody detects endogenous levels of total Aconitase 2
Immunogen Description	A synthesized peptide derived from human Aconitase 2
Other Names	ACO2; Aconitase2; Aconitate hydratase; ACONM; Citrate hydro lyase; ICRD;
Accession No.	Uniprot:Q99798
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:5000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50

Images



Western blot analysis of Aconitase 2 expression in K562 cell lysate.

Background

Catalyzes the isomerization of citrate to isocitrate via cis-aconitate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.