

Syntaxin 3 Rabbit mAb

Catalog No: #59899

Package Size: #59899-1 50ul #59899-2 100ul


Description

Product Name	Syntaxin 3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	Syntaxin 3 Antibody detects endogenous levels of total Syntaxin 3
Immunogen Description	A synthesized peptide derived from human Syntaxin 3
Other Names	Stx3; STX3A; Syntaxin 3A;
Accession No.	Uniprot:Q13277
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of Syntaxin 3 expression in HepG2 cell lysate.

Background

Potentially involved in docking of synaptic vesicles at presynaptic active zones.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.