

Fbx32 Rabbit mAb

Catalog No: #59909

Package Size: #59909-1 50ul #59909-2 100ul


Description

Product Name	Fbx32 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	Fbx32 Antibody detects endogenous levels of total Fbx32
Immunogen Description	A synthesized peptide derived from human Fbx32
Other Names	Atrogin 1; ATROGIN1; fbxo25; FBXO32; MAFbx;
Accession No.	Uniprot:Q969P5
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200

Images



Western blot analysis of Fbx32 expression in Human skeletal muscle lysate.

Background

Substrate recognition component of a (SKP1-CUL1-F-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.