

H Cadherin Rabbit mAb

Catalog No: #59913

Package Size: #59913-1 50ul #59913-2 100ul


Description

Product Name	H Cadherin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	H Cadherin Antibody detects endogenous levels of total H Cadherin
Immunogen Description	A synthesized peptide derived from human H Cadherin
Other Names	CDH13; CDHH; H-cadherin; Heart cadherin; P105; T cad; T Cadherin; Truncated cadherin;
Accession No.	Uniprot:P55290
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Images



Western blot analysis of H Cadherin expression in Human fetal heart lysate.

Background

Cadherins are calcium dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types. May act as a negative regulator of neural cell growth.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.