

GAB1 Rabbit mAb

Catalog No: #59919

Package Size: #59919-1 50ul #59919-2 100ul


Description

Product Name	GAB1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	GAB1 Antibody detects endogenous levels of total GAB1
Immunogen Description	A synthesized peptide derived from human GAB1
Other Names	GAB1; GRB 2 associated binder 1;
Accession No.	Uniprot:Q13480
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 FC 1:500

Images



Western blot analysis of GAB1 expression in HEK293 cell lysate.

Background

Probably involved in EGF and insulin receptor signaling.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.