

HIPK2 Rabbit mAb

Catalog No: #59951

Package Size: #59951-1 50ul #59951-2 100ul


Description

Product Name	HIPK2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human Mouse
Specificity	HIPK2 Antibody detects endogenous levels of total HIPK2
Immunogen Description	A synthesized peptide derived from human HIPK2
Other Names	hHIPk2; Hipk2; Nbak1; Stank;
Accession No.	Uniprot:Q9H2X6
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 IP 1:30 FC 1:500

Images



Western blot analysis of HIPK2 expression in HT-1080 cell lysate.

Background

Protein kinase acting as a corepressor of several transcription factors, including SMAD1 and POU4F1/Brn3a and probably NK homeodomain transcription factors. Inhibits cell growth and promotes apoptosis. Involved in transcriptional activation of TP53 and TP73.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.