

## ICAD Rabbit mAb

Catalog No: #59967

Package Size: #59967-1 50ul #59967-2 100ul


## Description

|                       |                                                                                                    |
|-----------------------|----------------------------------------------------------------------------------------------------|
| Product Name          | ICAD Rabbit mAb                                                                                    |
| Host Species          | Rabbit                                                                                             |
| Clonality             | Monoclonal                                                                                         |
| Isotype               | Rabbit IgG                                                                                         |
| Purification          | Affinity-chromatography                                                                            |
| Applications          | WB IHC ICC/IF IP FC                                                                                |
| Species Reactivity    | Human                                                                                              |
| Specificity           | ICAD Antibody detects endogenous levels of total ICAD                                              |
| Immunogen Description | A synthesized peptide derived from human ICAD                                                      |
| Other Names           | DFF1; DFF35; DFF45; DFFA; H13; ICAD; ICAD L; ICAD S;                                               |
| Accession No.         | Uniprot:O00273                                                                                     |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage               | Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.                     |

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50 FC 1:50

## Images



Western blot analysis of ICAD expression in HeLa cell lysate.

## Background

Inhibitor of the caspase-activated DNase (DFF40).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.