

PKA R2 Rabbit mAb

Catalog No: #59968

Package Size: #59968-1 50ul #59968-2 100ul


Description

Product Name	PKA R2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Cow
Specificity	PKA R2 Antibody detects endogenous levels of total PKA R2
Immunogen Description	A synthesized peptide derived from human PKA R2
Other Names	KAP2; PKR2; PRKAR2; PRKAR2A;
Accession No.	Uniprot:P13861
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50

Images



Western blot analysis of PKA R2 expression in K562 cell lysate.

Background

Regulatory subunit of the cAMP-dependent protein kinases involved in cAMP signaling in cells. Type II regulatory chains mediate membrane association by binding to anchoring proteins, including the MAP2 kinase.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.