

EWSR1 Rabbit mAb

Catalog No: #59969

Package Size: #59969-1 50ul #59969-2 100ul


Description

Product Name	EWSR1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC
Species Reactivity	Human Mouse Rat
Specificity	EWSR1 Antibody detects endogenous levels of total EWSR1
Immunogen Description	A synthesized peptide derived from human EWSR1
Other Names	EWS; Ewsr1; RNA binding protein EWS;
Accession No.	Uniprot:Q01844
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200

Images



Western blot analysis of EWSR1 expression in K562 cell lysate.

Background

Might normally function as a repressor. EWS-fusion-proteins (EFPS) may play a role in the tumorigenic process. They may disturb gene expression by mimicking, or interfering with the normal function of CTD-POLII within the transcription initiation complex.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.