

GSTM1 Rabbit mAb

Catalog No: #59976

Package Size: #59976-1 50ul #59976-2 100ul


Description

Product Name	GSTM1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Glutathione S Transferase mu Antibody detects endogenous levels of total Glutathione S Transferase mu
Immunogen Description	A synthesized peptide derived from human Glutathione S Transferase mu
Other Names	MU; H-B; GST1; GTH4; GTM1; MU-1; GSTM1-1; MGC26563; GSTM1a-1a; GSTM1b-1b;
Accession No.	Uniprot:P09488
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of GSTM1 expression in SH-SY5Y cell lysate.

Background

Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.