

SKP2 Rabbit mAb

Catalog No: #59982

Package Size: #59982-1 50ul #59982-2 100ul


Description

Product Name	SKP2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human
Specificity	SKP2 Antibody detects endogenous levels of total SKP2
Immunogen Description	A synthesized peptide derived from human SKP2
Other Names	FBL1; FBXL1; FLB1; p45; p45skp2; Skp2;
Accession No.	Uniprot:Q13309
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200

Images



Western blot analysis of SKP2 expression in HeLa cell lysate.

Background

Substrate recognition component of a SCF (SKP1-CUL1-F-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins involved in cell cycle progression, signal transduction and transcription.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.