

NDUFS1 Rabbit mAb

Catalog No: #60001

Package Size: #60001-1 50ul #60001-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NDUFS1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	NDUFS1 Antibody detects endogenous levels of total NDUFS1
Immunogen Description	A synthesized peptide derived from human NDUFS1
Other Names	CI-75kD; NDUFS1; PRO1304;
Accession No.	Uniprot:P28331
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:5000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50 FC 1:50

Images



Western blot analysis of NDUFS1 expression in Raji cell lysate.

Background

Core subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I) that is believed to belong to the minimal assembly required for catalysis. Complex I functions in the transfer of electrons from NADH to the respiratory chain.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.