

ITPA Rabbit mAb

Catalog No: #60078

Package Size: #60078-1 50ul #60078-2 100ul


Description

Product Name	ITPA Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	ITPA Antibody detects endogenous levels of total ITPA
Immunogen Description	A synthesized peptide derived from human ITPA
Other Names	Inosine triphosphatase; Itpa; ITPase; My049; NTPase;
Accession No.	Uniprot:Q9BY32
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Images



Western blot analysis of ITPA expression in HepG2 cell lysate.

Background

Hydrolyzes ITP and dITP to their respective monophosphate derivatives.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.