

MTAP Rabbit mAb

Catalog No: #60131

Package Size: #60131-1 50ul #60131-2 100ul

Description

Product Name	MTAP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	MTAP Antibody detects endogenous levels of total MTAP
Immunogen Description	A synthesized peptide derived from human MTAP
Other Names	BDMF; DMSFH; DMSMFH; LGMBF; MeSAdo phosphorylase; MSAP; MTAP; MTAPase;
Accession No.	Uniprot:Q13126
Calculated MW	31kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50 FC 1:50

Background

Plays a major role in polyamine metabolism and is important for the salvage of both adenine and methionine.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.