

## NgR3 Rabbit mAb

Catalog No: #60165

Package Size: #60165-1 50ul #60165-2 100ul


## Description

|                       |                                                                                                    |
|-----------------------|----------------------------------------------------------------------------------------------------|
| Product Name          | NgR3 Rabbit mAb                                                                                    |
| Host Species          | Rabbit                                                                                             |
| Clonality             | Monoclonal                                                                                         |
| Isotype               | Rabbit IgG                                                                                         |
| Purification          | Affinity-chromatography                                                                            |
| Applications          | WB IHC ICC/IF                                                                                      |
| Species Reactivity    | Human                                                                                              |
| Specificity           | NgR3 Antibody detects endogenous levels of total NgR3                                              |
| Immunogen Description | A synthesized peptide derived from human NgR3                                                      |
| Other Names           | NgR3; NGRH2; NGRL2; RTN4RL1;                                                                       |
| Accession No.         | Uniprot:Q86UN2                                                                                     |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage               | Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.                     |

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200

## Images



Western blot analysis of NgR3 expression in K562 cell lysate.

## Background

May play a role in regulating axonal regeneration and plasticity in the adult central nervous system.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.