

TEMT Rabbit mAb

Catalog No: #60176

Package Size: #60176-1 50ul #60176-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TEMT Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	TEMT Antibody detects endogenous levels of total TEMT
Immunogen Description	A synthesized peptide derived from human TEMT
Other Names	Amine N methyltransferase; Inmt; TEMT;
Accession No.	Uniprot:O95050
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of TEMT expression in A549 cell lysate.

Background

Catalyzes the N-methylation of tryptamine and structurally related compounds.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.