

KCTD9 Rabbit mAb

Catalog No: #60179

Package Size: #60179-1 50ul #60179-2 100ul


Description

Product Name	KCTD9 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human
Specificity	KCTD9 Antibody detects endogenous levels of total KCTD9
Immunogen Description	A synthesized peptide derived from human KCTD9
Other Names	BTBD27;
Accession No.	Uniprot:Q7L273
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200

Images



Western blot analysis of KCTD9 expression in HepG2 cell lysate.

Background

KCTD9 contains a potassium channel tetramerisation domain. The N-terminal, cytoplasmic tetramerisation domain (T1) of voltage-gated potassium channels encodes molecular determinants for subfamily-specific assembly of alpha-subunits into functional tetrameric channels.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.