

MDH2 Rabbit mAb

Catalog No: #60180

Package Size: #60180-1 50ul #60180-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MDH2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	MDH2 Antibody detects endogenous levels of total MDH2
Immunogen Description	A synthesized peptide derived from human MDH2
Other Names	M MDH; Malate dehydrogenase;
Accession No.	Uniprot:P40926
Calculated MW	36kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Background

Belongs to the LDH/MDH superfamily. MDH type 1 family.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.