

MLD Rabbit mAb

Catalog No: #60193

Package Size: #60193-1 50ul #60193-2 100ul

Description

Product Name	MLD Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	MLD Antibody detects endogenous levels of total MLD
Immunogen Description	A synthesized peptide derived from human MLD
Other Names	DEGS 1; DEGS; DES1; FADS7; MIG15; MLD;
Accession No.	Uniprot:O15121
Calculated MW	34kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Background

DEGS1 is a member of the membrane fatty acid desaturase family which is responsible for inserting double bonds into specific positions in fatty acids. It contains three His containing consensus motifs that are characteristic of a group of membrane fatty acid desaturases.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.