

MYO1B Rabbit mAb

Catalog No: #60201

Package Size: #60201-1 50ul #60201-2 100ul


Description

Product Name	MYO1B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	MYO1B Antibody detects endogenous levels of total MYO1B
Immunogen Description	A synthesized peptide derived from human MYO1B
Other Names	MMI alpha; MMIa; MYH 1c; Myo1b; myr1;
Accession No.	Uniprot:O43795
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Images



Western blot analysis of MYO1B expression in (1) 293 cell lysate; (2) Mouse heart lysate.

Background

Motor protein that may participate in process critical to neuronal development and function such as cell migration, neurite outgrowth and vesicular transport.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.