

ACAP2 Rabbit mAb

Catalog No: #60203

Package Size: #60203-1 50ul #60203-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ACAP2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	ACAP2 Antibody detects endogenous levels of total ACAP2
Immunogen Description	A synthesized peptide derived from human ACAP2
Other Names	ACAP2; Centaurin beta2; CENTB2; Cntb2;
Accession No.	Uniprot:Q15057
Calculated MW	88kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Background

GTPase-activating protein (GAP) for ADP ribosylation factor 6 (ARF6).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.