

CRCP Rabbit mAb

Catalog No: #60228

Package Size: #60228-1 50ul #60228-2 100ul


Description

Product Name	CRCP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	CRCP Antibody detects endogenous levels of total CRCP
Immunogen Description	A synthesized peptide derived from human CRCP
Other Names	CGRPRCP; CRCP; HsC17; RCP; RCP9;
Accession No.	Uniprot:O75575
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50

Images



Western blot analysis of CRCP expression in HeLa cell lysate.

Background

DNA-dependent RNA polymerase catalyzes the transcription of DNA into RNA using the four ribonucleoside triphosphates as substrates. Specific peripheric component of RNA polymerase III which synthesizes small RNAs, such as 5S rRNA and tRNAs.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.